

Multifit Pipette Tips



| Volume | Sterile | Packaging | Cat. No. |
|-----------|---------|---------------------------------|----------|
| 0.5-10 µl | No | Bulk | 310-2000 |
| 0.5-10 µl | No | Reload | 310-3000 |
| 0.5-10 µl | No | Racked | 310-4000 |
| 0.5-10 µl | Yes | Racked | 310-4050 |
| 0.5-10 µl | Yes | Low Retention, Racked | 310-4150 |
| 0.5-10 µl | No | Extended, Bulk | 320-2000 |
| 0.5-10 µl | No | Extended, Reload | 320-3000 |
| 0.5-10 µl | No | Extended, Racked | 320-4000 |
| 0.5-10 µl | Yes | Extended, Racked | 320-4050 |
| 0.5-10 µl | Yes | Extended, Low Retention, Racked | 320-4150 |
| 0.5-10 µl | No | Extended, Bulk | 325-2000 |
| 0.5-10 µl | No | Extended, Reload | 325-3000 |
| 0.5-10 µl | No | Extended, Racked | 325-4000 |
| 0.5-10 µl | Yes | Extended, Racked | 325-4050 |
| 0.5-10 µl | Yes | Extended, Low Retention, Racked | 325-4150 |
| 200 µl | No | Bulk | 340-2000 |
| 200 µl | No | Reload | 340-3000 |
| 200 µl | No | Racked | 340-4000 |
| 200 µl | Yes | Racked | 340-4050 |
| 200 µl | Yes | Low Retention, Racked | 340-4150 |
| 200 µl | No | Bulk | 345-2000 |
| 200 µl | No | Reload | 345-3000 |
| 200 µl | No | Racked | 345-4000 |
| 200 µl | Yes | Racked | 345-4050 |
| 200 µl | Yes | Low Retention, Racked | 345-4150 |
| 300 µl | No | Bulk | 350-2000 |
| 300 µl | No | Reload | 350-3000 |
| 300 µl | No | Racked | 350-4000 |
| 300 µl | Yes | Racked | 350-4050 |
| 300 µl | Yes | Low Retention, Racked | 350-4150 |
| 1250 µl | No | Extended, Bulk | 360-2000 |
| 1250 µl | No | Extended, Reload | 360-3000 |
| 1250 µl | Yes | Extended, Racked | 360-4000 |
| 1250 µl | Yes | Extended, Racked | 360-4050 |
| 1250 µl | Yes | Extended, Low Retention, Racked | 360-4150 |

10 µl, Short

10 µl, Extended

10 µl SlimFit, Extended

200 µl, PrecisePoint

200 µl Plus, PrecisePoint

300 µl, PrecisePoint

1250 µl, Extended

Packaging Options

Bulk

- Bottom gusset allows for a free-standing and stable bag
- Constructed from double-thick plastic
- Slider Seal provides effortless opening/ dosing

Reload System

- Designed to minimize waste and reduce ecological footprint
- Stacked inserts and tips are "reloaded" into empty racks
- Tip insert is "locked" into base for added security



Filter Barrier Tips



| Volume | Sterile | Packaging | Cat. No. |
|-----------|---------|---------------------------------|----------|
| 0.5-10 µl | Yes | Racked | 311-4050 |
| 0.5-10 µl | Yes | Low Retention, Racked | 311-4150 |
| 0.5-10 µl | No | Bulk | 311-2000 |
| 0.5-10 µl | Yes | Extended, Racked | 321-4050 |
| 0.5-10 µl | Yes | Extended, Low Retention, Racked | 321-4150 |
| 0.5-10 µl | No | Extended, Bulk | 321-2000 |
| 0.5-10 µl | Yes | Extended, Racked | 326-4050 |
| 0.5-10 µl | Yes | Extended, Low Retention, Racked | 326-4150 |
| 0.5-10 µl | No | Extended, Bulk | 326-2000 |
| 20 µl | Yes | Racked | 346-4050 |
| 20 µl | Yes | Low Retention, Racked | 346-4150 |
| 20 µl | No | Bulk | 346-2000 |
| 100 µl | Yes | Racked | 347-4050 |
| 100 µl | Yes | Low Retention, Racked | 347-4150 |
| 100 µl | No | Bulk | 347-2000 |
| 200 µl | Yes | Racked | 348-4050 |
| 200 µl | Yes | Low Retention, Racked | 348-4150 |
| 200 µl | No | Bulk | 348-2000 |
| 200 µl | Yes | Racked | 351-4050 |
| 200 µl | Yes | Low Retention, Racked | 351-4150 |
| 200 µl | No | Bulk | 351-2000 |
| 1000 µl | Yes | Extended, Racked | 361-4050 |
| 1000 µl | Yes | Extended, Low Retention, Racked | 361-4150 |
| 1000 µl | No | Extended, Bulk | 361-2000 |

10 µl, Short

10 µl, Extended

10 µl SlimFit, Extended

200 µl Plus, PrecisePoint

100 µl Plus, PrecisePoint

200 µl Plus, PrecisePoint

200 µl, PrecisePoint

1000 µl, Extended

Packaging Options

- Bulk:** 1000 tips/self-standing bag, 10 bags/case
- Racked:** 96 tips/rack, 10 racks/pack, 5 packs/case
- Reload:** 96 tips/tray, 10 trays/pack, 10 packs/case



E-Beam Sterilization

BioPointe pre-sterilized products are sterilized with an Electron Beam process. Unlike Gamma sterilization, E-Beam will not alter the color, or embrittle the plastic.



Racked

- Rigid construction ensures maximum stability with multi-channel pipetting
- 96-well, fast-snap insert features exact-hole alignment
- Snap-lock clasp keeps lid secure
- Recyclable, autoclavable, and individually numbered for product traceability

Pipette Tip Compatibility



| | Biohit ProLine / ProLine Plus / m-Line (Single) | Eppendorf Research Single, 8 and 12 channel | Eppendorf Reference Single channel | Finnpipette F2 Single channel | Pipetman Classic / Neo Single channel | HILDiscovery Comfort/Pro Labmate / Axyzet (Single) | Rainin Classic / XLS (standard shafts) | Nichiryo Single, 8&12 channel | Socorex Accura Single channel | VWR Signature / UHP Single channel |
|------------------------|---|---|------------------------------------|-------------------------------|---------------------------------------|--|--|-------------------------------|-------------------------------|------------------------------------|
| 0.1-2.5µl (PreLine) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.1-3µl (mPlus) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.5-10µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2-20µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 5-50µl (PreLine) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 10-100µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 20-200µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 100-1000µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.1-2.5µl (UltraMicro) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.5-10µl (UltraMicro) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 10-100µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 20-200µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 30-300µl (8/12ch) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 100-1000µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.1-2.5µl (UltraMicro) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.5-10µl (UltraMicro) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2-20µl (UltraMicro) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2-20µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 10-100µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 20-200µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 100-1000µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.2-2µl (Micro) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 1-10µl (Micro) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2-20µl (Micro) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2-20µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 10-100µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 20-200µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 100-1000µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.2-2µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 1-10µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2-20µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 20-100µl (Classic) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 50-200µl (Classic) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 200-1000µl (Classic) | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.2-2µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.5-10µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2-20µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 10-100µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 20-200µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 100-1000µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.1-2µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.5-10µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2-20µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 10-100µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 20-200µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 100-1000µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.1-2µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.5-10µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2-20µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 10-100µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 20-200µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 100-1000µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.1-2µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.5-10µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2-20µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 10-100µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 20-200µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 100-1000µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.1-2µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 0.5-10µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 2-20µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 10-100µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 20-200µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| 100-1000µl | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |

Multifit

10 µl, Short

10 µl, Extended

10 µl, SlimFit

200/200 µl Plus

300 µl

1250 µl

Filter Barrier

10 µl, Short

10 µl, Extended

10 µl, SlimFit

20 µl

100 µl

200 µl Plus¹

200 µl

1000 µl

For the full range of pipette tips and packaging options, please check our online catalog at www.biopointescientific.com.

This table is intended to be a general guide only. Pipette manufacturers occasionally change nose cone/shaft design and so tip fit can not be guaranteed.

¹ maximum volume with filter, 200µl

¹ maximum volume with filter, 200µl
©2018, all rights reserved, BioPointe Scientific Inc. PCR is a patented process of Roche Diagnostics. All brand or trade names are trademarks or registered trademarks of their respective holders.

Screw Cap Vials

| Volume | Style | Sterile | Cap Color | Packaging | Cat. No. |
|---------------------------|---------|---------|-----------|-----------|------------|
| 0.75 ml <i>Skirted</i> | Skirted | No | Natural | A | SBS75NA-BP |
| | Skirted | No | Rainbow | B | SBS75AS-BP |
| | Skirted | No | Color | A | SBS75**-BP |
| | Skirted | Yes | Natural | C | SBS75NA-PS |
| | Skirted | Yes | Rainbow | D | SBS75AS-PS |
| | Skirted | Yes | Color | C | SBS75**-PS |
| 1.6 ml <i>Skirted</i> | Skirted | No | Natural | A | SBS16NA-BP |
| | Skirted | No | Rainbow | B | SBS16AS-BP |
| | Skirted | No | Color | A | SBS16**-BP |
| | Skirted | Yes | Natural | C | SBS16NA-PS |
| | Skirted | Yes | Rainbow | D | SBS16AS-PS |
| | Skirted | Yes | Color | C | SBS16**-PS |
| 1.6 ml <i>Conical</i> | Conical | No | Natural | A | SBC16NA-BP |
| | Conical | No | Rainbow | B | SBC16AS-BP |
| | Conical | No | Color | A | SBC16**-BP |
| | Conical | Yes | Natural | C | SBC16NA-PS |
| | Conical | Yes | Rainbow | D | SBC16AS-PS |
| | Conical | Yes | Color | C | SBC16**-PS |
| 2.0 ml <i>Skirted</i> | Skirted | No | Natural | A | SBS20NA-BP |
| | Skirted | No | Rainbow | B | SBS20AS-BP |
| | Skirted | No | Color | A | SBS20**-BP |
| | Skirted | Yes | Natural | C | SBS20NA-PS |
| | Skirted | Yes | Rainbow | D | SBS20AS-PS |
| | Skirted | Yes | Color | C | SBS20**-PS |
| 2.0 ml <i>Conical</i> | Conical | No | Natural | A | SBC20NA-BP |
| | Conical | No | Rainbow | B | SBC20AS-BP |
| | Conical | No | Color | A | SBC20**-BP |
| | Conical | Yes | Natural | C | SBC20NA-PS |
| | Conical | Yes | Rainbow | D | SBC20AS-PS |
| | Conical | Yes | Color | C | SBC20**-PS |

- Perfect for long-term, secure storage of biological reagents
- Suitable for high speed centrifugation up to 17,000g
- Suitable for use at extreme temperatures ranging from -195°C (LN2 vapor phase) to +121°C

BioStar™ Screw-Cap Vials serve a wide range of life science laboratory needs, including sample collection, processing, centrifugation, long-term storage and packaging/shipping of biological reagents.

Select color of choice by using color code: **RD=Red, **YL**=Yellow, **GN**=Green, **BL**=Blue, **VT**=Violet, **AM**=Amber, **DA**=Dark Amber (For color selections AM and DA, both vials and caps are supplied in specified color)


RainbowPack = Assorted Cap Colors comprising Red, Yellow, Green, Blue and Violet, 2 Bags of each Cap Color

ClickSeal™ Microtubes

| Volume | Sterile | Type | Color | Packaging | Cat. No. |
|---------|---------|---------------|---------|-----------|-------------|
| 0.65 ml | No | Standard | Clear | A | SCN0645-BP |
| 1.7 ml | No | Standard | Clear | B | SCN1750-BP |
| 1.7 ml | No | Low Retention | Clear | C | SCN1750L-BP |
| 1.7 ml | No | Standard | Rainbow | D | SCA1750-BP |
| 2.0 ml | No | Standard | Clear | E | SCN2050-BP |

*Select color of choice by using color code:
R=Red, **Y**=Yellow, **G**=Green, **B**=Blue, **V**=Violet

MicroTube Rack System




Packaging Description
A = 960 Tubes/Bag, 5 Bags/Case
B = 480 Tubes/Bag, 5 Bags/Case
C = 120 Strips/Bag, 5 Bags/Case
D = 80 Strips/Bag, 5 Bags/Case
E = 96 Tubes/Rack, 10 Racks/Pack, 5 Packs/Case
F = 12 Strips/Rack, 10 Racks/Pack, 5 Packs/Case
G = 8 Strips/Rack, 10 Racks/Pack, 5 Packs/Case
H = 10 Mats/Bag, 5 Bags/Case
I = 96 Tubes/Rack, 6 Racks/Pack, 5 Packs/Case

PCR® Consumables


| Volume | Color | Description | Cat. No. |
|---------|-------|---|------------------|
| 0.2 ml | Clear | Domed | SPCST02-DC-01B |
| 0.2 ml | Clear | Flat Top | SPCST02-FF-01B |
| 0.2 ml | Clear | Flat Top, Optically Clear | SPCST02-RT-01B |
| 0.5 ml | Clear | Flat Top | SPCST05-FF-0 |
| 0.1 ml | Clear | 8-CapStrips™, Domed | SPCCS00-DC-08C |
| 0.2 ml | Clear | 8-CapStrips™, Flat Optically Clear | SPCCS00-RT-08C |
| 0.1 ml | Clear | 8-TubeStrips™, Low Profile | SPCTS01-LP-08B |
| 0.2 ml | Clear | 8-TubeStrips™ | SPCTS02-SP-08B |
| 0.2 ml | Clear | 8-TubeStrips™ and Domed CapStrips™ | SPCTS02-DC-08A |
| 0.2 ml | Clear | 8-TubeStrips™ and Optically Clear CapStrips™ | SPCTS02-RT-08A |
| 0.2 ml | Clear | 8-TubeStrips™ with Individually Attached Frosted Domed Caps | SPCTS02-IADC-08 |
| 0.2 ml | Clear | 8-TubeStrips™ with Individually Attached Flat Optically Clear RT Caps | SPCTS02-I AFC-08 |
| 0.1 ml | Clear | 8-TubeStrips™, Low Profile and Optically Clear CapStrips™ | SPCTS01-RT-08A |
| 0.65 ml | NS | A | STN0933-01B |
| 0.65 ml | NS | C | STN0933-08B |
| 0.65 ml | NS | E | STN0933-01R |
| 0.65 ml | S | E | STN0933-01RS |
| 0.65 ml | NS | F | STN0933-08R |
| 0.65 ml | S | F | STN0933-08RS |
| 1.2 ml | NS | A | STN0946-01B |
| 1.2 ml | NS | C | STN0946-08B |
| 1.2 ml | NS | D | STN0946-12B |
| 1.2 ml | NS | E | STN0946-01R |
| 1.2 ml | S | E | STN0946-01RS |
| 1.2 ml | NS | F | STN0946-08R |
| 1.2 ml | S | F | STN0946-08RS |
| 1.2 ml | NS | G | STN0946-12R |
| 1.2 ml | S | G | STN0946-12RS |

Packaging Description
8-CapStrips, Vented, Natural NS C STN3346-08C
8-CapStrips, Vented, Natural S C STN3346-08CS
12-CapStrips, Vented, Natural NS D STN3346-12C
12-CapStrips, Vented, Natural S D STN3346-12CS
8-CapStrips, Non-Vented, Natural NS C STN3346-09C
8-CapStrips, Non-Vented, Natural S C STN3346-09CS
MicroTube Rack System™ 96-Plug Mats NS H STN3346-96C
MicroTube Rack System™ 96-Plug Mats S H STN3346-96CS

Thin-Walled PCR Tubes



8-TubeStrips/CapStrips



For cyclor compatibility, plus sealing mats and film options, see the BioPointe Scientific online catalog at www.biopointescientific.com



96-well PCR Plates

| Volume | Color | Description | Cat. No. |
|--------|-------|---|----------------|
| 200µl | Clear | No-Skirt, Standard Profile | SPCP096-NS-SP |
| 100µl | Clear | No-Skirt, Low Profile | SPCP096-NS-LP* |
| 100µl | Clear | Half Skirt, Low Profile | SPCP096-HSR-LP |
| 200µl | Clear | Half Skirt, Standard Profile | SPCP096-HSR-SP |
| 100µl | Clear | Full Skirt, Low Profile | SPCP096-FS-LP |
| 200µl | Clear | Half Skirt, Standard Profile | SPCP096-HS-SP |
| 100µl | Clear | Full Skirt, for Roche® LightCycler® 480 | SPCP096-L480 |
| 100µl | White | Full Skirt, for Roche® LightCycler® 480 | SPWP096-L480 |

384-well PCR Plates

| | | | |
|------|-------|---|---------------|
| 50µl | Clear | Full Skirt, Single Notch | SPCP384-FS-1N |
| 50µl | Clear | Full Skirt, Dual Notch | SPCP384-FS-2N |
| 50µl | Clear | Full Skirt, for Roche® LightCycler® 480 | SPCP384-L480 |
| 50µl | White | Full Skirt, for Roche® LightCycler® 480 | SPWP384-L480 |

Packaging Description

Tubes: 1000/Pack, 10 Packs/Case
Strips: 125 Strips/Pack, 10 Packs/Case

Plates: 10 Plates/Pack, 10 Packs/Case
 *20 Plates/Pack, 5 Packs/Case



FOR A COMPLETE ONLINE CATALOG OF PRODUCTS, VISIT OUR WEBSITE AT WWW.BIOPINTESCIENTIFIC.COM